




VPS37A Antibody

Product Code	CSB-PA256766
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEZ2
Immunogen	Fusion protein of Human VPS37A
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	VPS37A


Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA256766(VPS37A Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using CSB-PA256766(VPS37A Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-2:
Human bladder carcinoma and fetal liver tissue,
Primary antibody: CSB-PA256766(VPS37A
Antibody) at dilution 1/300 dilution, Secondary
antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 15 seconds