




6xHis Tag Antibody

Product Code	CSB-PA270215
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Synthetic peptide: 18xHis conjugated to KLH
Raised In	Rabbit
Species Reactivity	HIS Tag
Tested Applications	ELISA,WB;ELISA:1:5000-1:10000,WB:1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG

Image



Gel: 8%SDS-PAGE, Lysate: 0.5 µg, Lane: Recombinant 6HIS Tagged protein , Primary antibody: CSB-PA270215(6His Antibody) at dilution 1/2000 (lane 1),1/5000 (lane 2),1/8000 (lane 3),1/10000 (lane 4), Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.