




# HA Tag Antibody

<b>Product Code</b>	CSB-PA275079
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Synthetic peptide: YPYDVDPDYA conjugated to KLH.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	HA Tag
<b>Tested Applications</b>	ELISA,WB;ELISA:1:2000-1:10000,WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG

## Image



Gel: 8%SDS-PAGE, Lysate: 2 µg, Lane: Recombinant HA tagged protein, Primary antibody: CSB-PA275079(HA Tag Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 second

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.