





ALK Antibody

Product Code	CSB-PA847171
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM73
Immunogen	Synthetic peptide of Human ALK
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,IHC;ELISA:1:1000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	ALK

Image



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA847171(ALK Antibody) at dilution 1/25. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA847171(ALK Antibody) at dilution 1/25. (Original magnification: ×200)