



NPPA Antibody

Product Code	CSB-PA578654
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01160
Immunogen	Synthetic peptide of Human NPPA
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA; ELISA: 1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	NPPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.