




DPYSL5 Antibody

Product Code	CSB-PA996118
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPU6
Immunogen	Synthetic peptide of Human DPYSL5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:500-1:5000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	DPYSL5

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane:
 Human fetal brain tissue, Primary antibody:
 CSB-PA996118(DPYSL5 Antibody) at dilution
 1/150, Secondary antibody: Goat anti rabbit IgG
 at 1/8000 dilution, Exposure time: 15 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.