




CCNB1 Antibody

Product Code	CSB-PA018543
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14635
Immunogen	Synthetic peptide of Human CCNB1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	CCNB1

Image



Gel: 8%SDS-PAGE,Lysate: 40 µg,Lane 1-3:
Jurkat cells, 231 cells, Hela cells,Primary
antibody: CSB-PA018543(CCNB1 Antibody) at
dilution 1/300 dilution,Secondary antibody: Goat
anti rabbit IgG at 1/8000 dilution,Exposure time:
10 seconds