




EIF4G2 Antibody

Product Code	CSB-PA063947
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78344
Immunogen	Synthetic peptide of Human EIF4G2
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	EIF4G2

Image



Gel: 6%SDS-PAGE,Lysate: 40 µg,Lane 1-2:
Hela cells, HepG2 cells,Primary antibody: CSB-
PA063947(EIF4G2 Antibody) at dilution 1/500
dilution,Secondary antibody: Goat anti rabbit IgG
at 1/8000 dilution,Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.