




PVRL1 Antibody

Product Code	CSB-PA139246
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15223
Immunogen	Synthetic peptide of Human PVRL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:500-1:5000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	NECTIN1

Image



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-3: Human cervical cancer tissue, Human legs fibrous histiocytoma, Human fetal brain tissue, Primary antibody: CSB-PA139246(PVRL1 Antibody) at dilution 1/400, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.