




FOXL1 Antibody

Product Code	CSB-PA236731
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12952
Immunogen	Synthetic peptide of Human FOXL1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	FOXL1

Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-2:
293T cells, mouse intestine tissue,
Primary antibody: CSB-PA236731(FOXL1
Antibody) at dilution 1/850, Secondary antibody:
Goat anti rabbit IgG at 1/8000 dilution, Exposure
time: 30 seconds