




HMGB2 Antibody

Product Code	CSB-PA750763
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26583
Immunogen	Synthetic peptide of Human HMGB2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:10000,WB:1:2000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	HMGB2

Image



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-2: Hela cells, 231 cells, Primary antibody: CSB-PA750763(HMGB2 Antibody) at dilution 1/2000, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 50 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.