





HMGB3 Antibody

Product Code	CSB-PA780376
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15347
Immunogen	Synthetic peptide of Human HMGB3
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	HMGB3

Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA780376(HMGB3 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-2: Mouse lung tissue, Mouse brain tissue, Primary antibody: CSB-PA780376(HMGB3 Antibody) at dilution 1/700, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1.5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.