




ING1 Antibody

Product Code	CSB-PA824082
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK53
Immunogen	Synthetic peptide of Human ING1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:500-1:1000,WB:1:100-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	ING1

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-4: Mouse brain tissue, Jurkat cells, Raw264.7 cells, hepG2 cells, Primary antibody: CSB-PA824082(ING1 Antibody) at dilution 1/53, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 second

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.