




C1QTNF5 Antibody

Product Code	CSB-PA031979
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXJ0
Immunogen	Synthetic peptide of Human C1QTNF5
Raised In	Guinea pig
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; ELISA: 1:1000-1:2000, WB: 1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	C1QTNF5

Image



Gel: 10+12%SDS-PAGE, Lysate: 40 µg, Lane: A549 cells, Primary antibody: CSB-PA031979(C1QTNF5 Antibody) at dilution 1/300, Secondary antibody: Donkey anti Guinea pig IgG at 1/8000 dilution, Exposure time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.