





RAD52 Antibody

Product Code	CSB-PA955701
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43351
Immunogen	Synthetic peptide of Human RAD52
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:2000-1:10000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	RAD52 homolog (S. cerevisiae)
Immunogen Species	Homo sapiens (Human)
Target Names	RAD52

Image



The image is immunohistochemistry of paraffin-embedded Human colon cancer using CSB-PA955701(RAD52 Antibody) at dilution 1/50. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer using CSB-PA955701(RAD52 Antibody) at dilution 1/50. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.