





SPAG4 Antibody

Product Code	CSB-PA967102
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPE6
Immunogen	Synthetic peptide of Human SPAG4
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:3000,WB:1:100-1:500,IHC:1:15-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	SPAG4

Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA967102(SPAG4 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 50 µg, Lane: NIH/3T3 cells, Primary antibody: CSB-PA967102(SPAG4 Antibody) at dilution 1/150, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 30 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.