





PRKD3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA145264 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O94806 |
| Immunogen | Synthetic peptide of Human PRKD3 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Tested Applications | ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PRKD3 |


Image



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA145264(PRKD3 Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA145264(PRKD3 Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: x200)



Gel: 6%SDS-PAGE, Lysate: 40 μ g, Lane 1-6:
 Human testis tissue, K562 cells, A549 cells, Raji
 cells, NIH/3T3 cells, Hela cells, Primary
 antibody: CSB-PA145264(PRKD3 Antibody) at
 dilution 1/200, Secondary antibody: Goat anti
 rabbit IgG at 1/8000 dilution, Exposure time: 10
 minutes