




PAK3 Antibody

Product Code	CSB-PA103384
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75914
Immunogen	Synthetic peptide of Human PAK3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	PAK3

Image



Gel: 6%SDS-PAGE,Lysate: 40 µg,Lane 1-2:
Human fetal brain tissue, Mouse brain
tissue,Primary antibody: CSB-PA103384(PAK3
Antibody) at dilution 1/200 dilution,Secondary
antibody: Goat anti rabbit IgG at 1/8000
dilution,Exposure time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.