




ARTN Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA167947 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5T4W7 |
| Immunogen | Synthetic peptide of Human ARTN |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ARTN |

Image



Gel: 12%SDS-PAGE, Lysate: 50 µg, Lane 1-3: 293T cells, SP20 cells, human fetal liver tissue, Primary antibody: CSB-PA167947(ARTN Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.