



Exendin-4 Antibody

Product Code	CSB-PA853477
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26349
Immunogen	Synthetic peptide of Heloderma suspectum Exendin 4
Raised In	Rabbit
Species Reactivity	Heloderma suspectum
Tested Applications	ELISA;ELISA:1:1000-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Heloderma suspectum (Gila monster)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.