





HDLBP Antibody

Product Code	CSB-PA876442
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00341
Immunogen	Synthetic peptide of Human HDLBP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	HDLBP


Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA876442(HDLBP Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: X200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA876442(HDLBP Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: X200)



Gel: 6%SDS-PAGE, Lysate: 40μg, Lane:
NIH/3T3 cells, Primary antibody: CSB-
PA876442(HDLBP Antibody) at dilution 1/750,
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution, Exposure time: 1 minute