





# AGPAT9 Antibody

<b>Product Code</b>	CSB-PA798773
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q53EU6?
<b>Immunogen</b>	Synthetic peptide of Human AGPAT9
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Tested Applications</b>	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:25-1:100
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GPAT3


## Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA798773(AGPAT9 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA798773(AGPAT9 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 80  $\mu$ g, Lane:  
Humna placenta tissue, Primary antibody: CSB-  
PA798773(AGPAT9 Antibody) at dilution 1/200  
dilution, Secondary antibody: Goat anti rabbit  
IgG at 1/8000 dilution, Exposure time: 15  
seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.