




NPHS2 Antibody

Product Code	CSB-PA579527
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP85
Immunogen	Synthetic peptide of Human NPHS2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000,
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification Method	Antigen affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NPHS2

Image



Western blot analysis of lysates from Jurkat cells, primary antibody was diluted at 1:1000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.