



PPARGC1A Antibody

Product Code	CSB-PA042134
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBK2
Immunogen	Synthetic peptide of Human PPARGC1A
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA; ELISA:1:1000-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPARGC1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.