



TNFRSF11A Antibody

Product Code	CSB-PA925190
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6Q6
Immunogen	Synthetic peptide of Human TNFRSF11A
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF11A

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-2:
Human liver tissue, mouse kidney tissue,
Primary antibody: CSB-PA925190(TNFRSF11A
Antibody) at dilution 1/200, Secondary antibody:
Goat anti rabbit IgG at 1/8000 dilution, Exposure
time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.