





HIF1AN Antibody

Product Code	CSB-PA959614
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWT6
Immunogen	Synthetic peptide of Human HIF1AN
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:10000,WB:1:500-1:2000,IHC:1:50-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	HIF1AN


Image



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using CSB-PA959614(HIF1AN Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA959614(HIF1AN Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



Gel: 8%SDS-PAGE, Lysate: 60 μ g, Lane 1-2:
 Mouse brain tissue, Mouse liver tissue, Primary
 antibody: CSB-PA959614(HIF1AN Antibody) at
 dilution 1/200, Secondary antibody: Goat anti
 rabbit IgG at 1/8000 dilution, Exposure time: 20
 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.