




Ly6a Antibody

Product Code	CSB-PA983699
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05533
Immunogen	Synthetic peptide of Mouse Ly6a
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Mus musculus (Mouse)
Target Names	Ly6a

Image



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane: Mouse skin tissue, Primary antibody: CSB-PA983699(Ly6a Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes