




PYY Antibody

Product Code	CSB-PA831475
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10082
Immunogen	Synthetic peptide of Human PYY
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	Antigen affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PYY

Image



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane:
 Mouse small intestines tissue, Primary antibody:
 CSB-PA831475(PYY Antibody) at dilution 1/200,
 Secondary antibody: Goat anti rabbit IgG at
 1/8000 dilution, Exposure time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.