




SPON1 Antibody

Product Code	CSB-PA546167
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCB6
Immunogen	Synthetic peptide of Human SPON1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:5000-1:20000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification Method	Antigen affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPON1

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane 1-2: Mouse brain tissue, Mouse lung tissue, Primary antibody: CSB-PA546167(SPON1 Antibody) at dilution 1/400, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.