




KCNK3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA964419 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O35111 |
| Immunogen | Synthetic peptide of Human KCNK3 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Mus musculus (Mouse) |
| Target Names | Kcnk3 |

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: Mouse heart tissue, Primary antibody: CSB-PA964419(KCNK3 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.