




TPD52L1 Antibody

Product Code	CSB-PA215944
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16890
Immunogen	Synthetic peptide of Human TPD52L1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	TPD52L1

Image



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane 1-5: 293T cells, 231 cells, PC3 cells, A549 cells, K562 cells, Primary antibody: CSB-PA215944(TPD52L1 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.