





VPS26A Antibody

Product Code	CSB-PA234289
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75436
Immunogen	Synthetic peptide of Human VPS26A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:30-1:150
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	VPS26A


Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA234289(VPS26A Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA234289(VPS26A Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



Gel: 12%SDS-PAGE,Lysate: 40 µg,Lane 1-7: 293T cells, A431 cells, Hela cells, HepG2 cells, Human fetal liver tissue, PC3 cells, Human placenta tissue,Primary antibody: CSB-PA234289(VPS26A Antibody) at dilution 1/200 dilution,Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,Exposure time: 40 seconds



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.