



E.coli Antibody

Product Code	CSB-PA152181
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Escherichia coli
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA;ELISA:1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Escherichia coli
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.