



SPAM1 Antibody, HRP conjugated

Product Code	CSB-PA022471LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38567
Immunogen	Recombinant Human Hyaluronidase PH-20 protein (36-490AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronidase PH-20 (Hyal-PH20) (EC 3.2.1.35) (Hyaluronoglucosaminidase PH-20) (Sperm adhesion molecule 1) (Sperm surface protein PH-20), SPAM1, HYAL3 PH20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SPAM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.