



LPA Antibody, FITC conjugated

Product Code	CSB-PA362446YC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08519
Immunogen	Recombinant Human Apolipoprotein(a) protein (3999-4112AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein(a) (Apo(a)) (Lp(a)) (EC 3.4.21.-), LPA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.