




# EIF1AX Antibody

<b>Product Code</b>	CSB-PA007506DSR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P47813
<b>Immunogen</b>	Recombinant Human Eukaryotic translation initiation factor 1A, X-chromosomal protein (10-144AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Eukaryotic translation initiation factor 1A, X-chromosomal (eIF-1A X isoform) (Eukaryotic translation initiation factor 4C) (eIF-4C), EIF1AX, EIF1A EIF4C
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	EIF1AX

## Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA007506DSR2HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.