





CD34 Antibody

Product Code	CSB-PA001451
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28906
Immunogen	Synthesized peptide derived from the Internal region of Human CD34.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD34; Hematopoietic progenitor cell antigen CD34; CD antigen CD34
Immunogen Species	Homo sapiens (Human)
Target Names	CD34


Image



Western Blot analysis of HepG2 cells using CD34 Polyclonal Antibody



Western Blot analysis of K562 HepG2 cells using CD34 Polyclonal Antibody



mouse-lung mouse-brain mouse-spleen mouse-kidney

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.