



SLC34A2 Antibody, Biotin conjugated

Product Code	CSB-PA021581LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95436
Immunogen	Recombinant Human Sodium-dependent phosphate transport protein 2B protein (574-689AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent phosphate transport protein 2B (Sodium-phosphate transport protein 2B) (Na(+)-dependent phosphate cotransporter 2B) (NaPi3b) (Sodium/phosphate cotransporter 2B) (Na+)/Pi cotransporter 2B) (NaPi-2b) (Solute carrier family 34 member 2), SLC34A2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC34A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.