



SLC34A2 Antibody, Biotin conjugated

Product Code	CSB-PA021581LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95436
Immunogen	Recombinant Human Sodium-dependent phosphate transport protein 2B protein (574-689AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent phosphate transport protein 2B (Sodium-phosphate transport protein 2B) (Na(+)-dependent phosphate cotransporter 2B) (NaPi3b) (Sodium/phosphate cotransporter 2B) (Na+)/Pi cotransporter 2B) (NaPi-2b) (Solute carrier family 34 member 2), SLC34A2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC34A2