



Recombinant Mycobacterium tuberculosis

Diacylglycerol acyltransferase/mycolyltransferase

Ag85C (fbpC)

Product Code	CSB-EP363614MVZ
Relevance	The antigen 85 proteins (FbpA, FbpB, FbpC) are responsible for the high affinity of mycobacteria to fibronectin, a large adhesive glycoprotein, which facilitates the attachment of M.tuberculosis to murine alveolar macrophages (AMs). They also help to maintain the integrity of the cell wall by catalyzing the transfer of mycolic acids to cell wall arabinogalactan and through the synthesis of alpha,alpha-trehalose dimycolate (TDM, cord factor). They catalyze the transfer of a mycoloyl residue from one molecule of alpha,alpha-trehalose monomycolate (TMM) to another TMM, leading to the formation of TDM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WQN8
Alias	Acyl-CoA:diacylglycerol acyltransferase Antigen 85 complex C Short name: 85C Short name: Ag85C Fibronectin-binding protein C Short name: Fbps C
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	AFSRPGLPVEYLQVPSASMGRDIKVQFQGGGPHAVYLLDGLRAQDDYNGWDI NTPAFEEYYQSGLSVIMPVGGQSSFYTDWYQPSQSNGQNYTYKWETFLTRE MPAWLQANKGVSPTGNAAVGLSMGGSALILAAAYYPQQFPYAASLSGFLNPS EGWWPTLIGLAMNDSGGYNANSMWGPSSDPAWKRNDPMVQIPRLVANNTRI WVYCGNGTPSDLGGDNIPAKFLEGLTLRTNQTFRDYAADGGRNGVFNFPPN GTHSWPYWNEQLVAMKADIQHVLNGATPPAAPAAPAA
Lead Time	3-7 business days
Research Area	Others
Source	E.coli
Gene Names	fbpC
Expression Region	46-340aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	36.1kDa

