




# Recombinant Mesomycoplasma hyorhinis High affinity transport system protein p37 (p37)

<b>Product Code</b>	CSB-EP323355MLR
<b>Relevance</b>	P37 is part of a high-affinity transport system.
<b>Abbreviation</b>	Recombinant Mesomycoplasma hyorhinis p37 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P15363
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mesomycoplasma hyorhinis (Mycoplasma hyorhinis)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	CSNTGVVKQEDVSVSQGQWDKSITFGVSEAWLNKKKGGEKVNKEVINTFLEN FKKEFNKLNANDKTKNFDDVDFKVTPIQDFTVLLNNLSTDNPELDFGINASGK LVEFLKNNPGIITPALETTTNSFVFDKEKDKFYVDGTDSDPLVKIAKEINKIFVET PYASWTDENHKWNGNVYQSVYDPTVQANFYRGMWIKGNDETLAKIKKAWN DKDWNTFRNFGILHGKDNSFSKFKLEETILKNHFQNKFTTLNEDRSAHPNAYK QKSADTLGTLDDFHIAFSEEGSFAWTHNKSATKPFETKANEMEALIVTNPIPY DVGVFRKSVNQLEQNLIVQTFINLAKNKQDITYGPLLGYNGYKKIDNFQKEIVEV YEKAIAK
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	p37
<b>Expression Region</b>	24-403aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	59.5kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.