





Recombinant Human Odorant-binding protein 2a (OBP2A)

Product Code	CSB-EP873688HU
Relevance	Probably binds and transports small hydrophobic volatile molecules with a higher affinity for aldehydes and large fatty acids.
Abbreviation	Recombinant Human OBP2A protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NY56
Alias	Odorant-binding protein IIa Short name:OBPIIa
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	LSFTLEEEDITGTWYVKAMVVDKDFPEDRRPRKVSPVKVTALGGGNLEATFTF MREDRCIQKKILMRKTEEPGKFSAYGGRKLIYLQELPGTDDYVFYCKDQRRGG LRYMGKLVGRNPNTNLEALEEFKKLVQHKGLSEEDIFMPLQTGSCVLEH
Research Area	Neuroscience
Source	E.coli
Target Names	OBP2A
Expression Region	16-170aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	33.8kDa
Protein Length	Full Length of Mature Protein
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.