




# Recombinant Pongo pygmaeus Phospholipid hydroperoxide glutathione peroxidase (GPX4) (U73S), partial

<b>Product Code</b>	CSB-EP670014EXP
<b>Relevance</b>	Could play a major role in protecting mammals from the toxicity of ingested lipid hydroperoxides. Essential for embryonic development. Protects from radiation and oxidative damage
<b>Abbreviation</b>	Recombinant ongo pygmaeus GPX4 protein (U73S), partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4AEH2
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	ongo pygmaeus (Bornean orangutan)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MSLGRLCRLLPALLCGALAAPGLAGTMCASRDDWRCARSMHEFSAKDIDGH MVNLDKYRGFVCIVTNVASQUQKTEVNYTQLVDLHARYAECGLRILAFPCNQF GKQEPGSNEEIKEFAAGYNVKFDMFSKICVNGDDAHPLWKWMKIQPKGKIL GNAIKWNFTKFLI
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	GPX4
<b>Expression Region</b>	28-197aa(U73S)
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	35.0kDa
<b>Protein Length</b>	Cytoplasmic Domain
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.