




Recombinant Human Pancreatic adenocarcinoma up-regulated factor (ZG16B)

Product Code	CSB-EP836195HU
Abbreviation	Recombinant Human ZG16B protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96DA0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GKMYGPGGGKYFSTTEDYDHEITGLRVSVGLLLVKSQVVKLGDSWDVKLGAL GGNTQEVTLQPGEYITKVFVAFQAFRLRGMVMYTSKDRYFYFGKLDGQISSAYP SQEGQVLVGIYGQYQLLGIKISIGFEWNYPLEEPTTEPPVNLTYSANSPVGR
Research Area	Immunology
Source	E.coli
Target Names	ZG16B
Expression Region	53-208aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	21.2kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.