




# Recombinant Mouse Complement receptor type 2 (Cr2), partial

<b>Product Code</b>	CSB-EP005934MO1
<b>Relevance</b>	Receptor for complement C3d. Participates in B lymphocytes activation.
<b>Abbreviation</b>	Recombinant Mouse Cr2 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P19070
<b>Alias</b>	Complement C3d receptor CD_antigen: CD21
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	ISCDPPPEVKNARKPYYSPLIVPGTVLRYTCSPSYRLIGEKAIFCISENQVHATW DKAPPICESVNKTISCSDPIVPGGFMNKGSKAPFRHGDSVTFTCKANFTMKGS KTVWCQANEMWGPTALPVCES
<b>Research Area</b>	Immunology
<b>Source</b>	E.coli
<b>Target Names</b>	Cr2
<b>Expression Region</b>	12-140aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	18.1kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.