




# Recombinant Mouse Beta-adrenergic receptor kinase 1 (Grk2)

<b>Product Code</b>	CSB-EP858153MO
<b>Relevance</b>	Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors, probably inducing a desensitization of them.
<b>Abbreviation</b>	Recombinant Mouse Grk2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99MK8
<b>Alias</b>	G-protein-coupled receptor kinase 2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MADLEAVLADVSYLMAMEKSKATPAARASKKILLPEPSIRSVMQKYLEDRGEV TFEKIFSQKLGYLLFRDFCLNHLEEAKPLVEFYEEIKKYEKLETEEERVVRSREI FDSYIMKELLACSHPFKSNATEHVQGHVKKQVPPDLFQPYIEEICQNLRGDVF QKFIESDKFTRFCQWKNVELNIHLTMNDFSVHRIIGRGGFGEVYGCRKADTGK MYAMKCLDKKRIKMKQGETLALNERIMLSLVSTGDPCFIVCMSYAFHTPKLSF ILDLMNGGDLHYHLSQHGVFSEADMRFYAAEIIILGLEHMHNRFVYRDLKPANI LLEDEHGHVRISDLGLACDFSKKRPHASVGTGTHGYMAPEVLQKGVAYDSSADWF SLGCMLFKLLRGHSPFRQHKTKDKHEIDRMTLTMAVELPDSFSPELRSLLEGL LQRDVNRRLGCLGRGAQEVKESPFRRSLDWQMVFLQKYPPLIPPRGEVNAA DAFDIGSFDEEDTKGIKLLDSDQELYRNFPLTISERWQQEVAETVFDTINAETD RLEARKKAKNKQLGHEEDYALGKDCIVHGYSKMGNPFLTQWQRRYFYLFPN RLEWRGEGEAPQSLLTMEEIQSVEETQIKERKCLLLKIRGGKQFVLQCDSDPE LVQWKELRDAYREAQQLVQRVPMKNKPRSPVVVELSKVPLIQRGSANGL
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	Grk2
<b>Expression Region</b>	1-689aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	95.6kDa
<b>Protein Length</b>	Full Length



## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.