






Recombinant Staphylococcus aureus Gamma-hemolysin component C (hlgC)

Product Code	CSB-EP739984SKX
Relevance	Toxin that seems to act by forming pores in the membrane of the cell. Has a hemolytic and a leucotoxic activity
Abbreviation	Recombinant Staphylococcus aureus hlgC protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GE13
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ANDTEDIGKGNDVEIIKRTEDKTSNKWGVQTQNIQFDFVKDKKYNKDALILKMQG FISSRTTYNYKNTNHIKSMRWPFQYNIGLKTNDKYVSLINYLPKNKIESTNVSQ TLGYNIGGNFQSAPSLGGNGSFNYSKSIYSTQQNYVSEVEQQNSKSVLWGVK ANSFATESGQKSAFDSDLFVGYPKPHSKDPRDYFVPDSELPPLVQSGFNPSFIA TVSHEKGSSDTSEFEITYGRNMDVTHAIKRSTHYGNSYLDGHRVHNAFKNRN YTVKYEVDNWKTHEIKVKGQN
Research Area	Others
Source	E.coli
Target Names	hlgC
Expression Region	30-315aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	48.6kDa
Protein Length	Full Length of Mature Protein
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.