




AAR2 Antibody

Product Code	CSB-PA896704LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y312
Immunogen	Recombinant Human Protein AAR2 homolog protein (2-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein AAR2 homolog (AAR2 splicing factor homolog), AAR2, C20orf4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AAR2

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA896704LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.