



# Abhd11 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA811717LA01MO  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8K4F5  |
| <b>Immunogen</b>           | Recombinant Mouse Protein ABHD11 protein (1-307AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Mouse   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Protein ABHD11 (EC 3.-.-.) (Alpha/beta hydrolase domain-containing protein 11) (Abhydrolase domain-containing protein 11) (Williams-Beuren syndrome chromosomal region 21 protein homolog), Abhd11, Wbscr21 |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)  |
| <b>Research Area</b>       | Others  |
| <b>Target Names</b>        | Abhd11  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |