




ACTR8 Antibody

Product Code	CSB-PA872523LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H981
Immunogen	Recombinant Human Actin-related protein 8 protein (35-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 8 (hArp8) (INO80 complex subunit N), ACTR8, ARP8 INO80N
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ACTR8

Image



Western blot
 All lanes: ACTR8 antibody at 16µg/ml
 Lane 1: A549 whole cell lysate
 Lane 2: HCT116 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 71, 59, 37 kDa
 Observed band size: 71 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.